

Version with Markings to Show Changes Made

The amended paragraphs indicate deletions by ~~strikeout~~ and insertions by underlining.

In the specification, page 3, second full paragraph.

Ballabio et al. (1991), disclose a ~~single-tube-multiplex~~
~~allele-specific~~ single-tube, multiplex allele-specific PCR test using
two different dye-tagged fluorescent primers for detection of the
ΔF508 cystic fibrosis mutation.

In the specification, page 11, first full paragraph.

The primers must also be designed so that the size of the
resulting amplification products differ in length, thereby facilitating
assignment of alleles to individual loci during detection.
Inappropriate selection of primers can produce several undesirable
effects such as lack of amplification, amplification at multiple sites,
primer dimer formation, undesirable interaction of primer
sequences from different loci, production of alleles from one locus
which overlap with alleles from another, or ~~requirement~~ the need
for amplification conditions or protocols for the different loci
which are incompatible in a multiplex. The synthesis of the
primers is conducted by procedures known to those skilled in the
art.

In the specification, page 18, third full paragraph.

In this example, a DNA template was amplified at the
individual loci HUMCSF1PO, HUMTPOX, HUMTH01, and
HUMVWFA31 simultaneously in a single reaction vessel. The
PCR amplifications were performed in 25μl volumes using 25ng
template, 0.04U *Taq* DNA Polymerase/μl, 1x STR Buffer (50mM
KCl, 10mM Tris-HCl (pH 9.0 at 25°C), 0.1% Triton X-100, 1.5mM
MgCl₂ and 200μM each of dATP, dCTP, dGTP and dTTP), and
using a Thermal Cycler 480 (Perkin Elmer Cetus). Amplification
protocol 1, as described in Example 1, was employed. Eight

amplification primers were used in combination, including 1 μ M each HUMCSF1PO primer 2 [~~SEQ. ID. 5~~] [SEQ. ID. 6] and fluorescein-labeled primer 1 [SEQ. ID. 5], 0.15 μ M each HUMTPOX primer 1 [SEQ. ID. 29] and fluorescein-labeled primer 2 [SEQ. ID. 30], 0.2 μ M each HUMTH01 primer 2 [SEQ. ID. 28] and fluorescein-labeled primer 1 [SEQ. ID. 27], and 1 μ M each HUMVWFA31 primer 1 [SEQ. ID. 31] and fluorescein-labeled primer 2 [SEQ. ID. 32].